



PUR-Blue™ Swab Sampler

A versatile tool for better sample collection



Unlike traditional fiber tipped swabs, PUR-Blue™ Swab Samplers use medical grade polyurethane foam tips, that resist tearing and fraying, making it a better choice for sampling rough surfaces. Manufactured to consistent rigorous standards, polyurethane is non-toxic and biocide free, which greatly improves the recovery and release of collected microorganisms, without any components that may interfere with downstream test methods. Try it with HiCap™ (High Capacity) Neutralizing Broth, the first collection broth developed specifically for use in the food and beverage industry. By neutralizing residual sanitizers on the sampled surface and maintaining the viability of the collected microorganisms for up to 72 hours, HiCap prevents false negative test results and accurately quantifies microorganisms.



Get Better Results

For More Information
Call 877-260-6441 or Email
info@worldbioproducts.net

Large Tip for Qualitative Testing, 100/case	Item #
HiCap™ Neutralizing Broth, 1.5 ml	BLU-HC-P
D/E Neutralizing Broth, 1.5 ml	BLU-DE-P
Lethen Broth, 1.5 ml	BLU-LET-P
Neutralizing Buffer, 1.5 ml	BLU-NB-P

Large Tip for Quantitative Testing, 100 case	
HiCap™ Neutralizing Broth, 5 ml	BLU-5HC
HiCap™ Neutralizing Broth, 10 ml	BLU-10HC
D/E Neutralizing Broth, 10 ml	BLU-10DE **
Lethen Broth, 5 ml	BLU-5LET
Lethen Broth, 10 ml	BLU-10LET
Neutralizing Buffer, 4 ml	BLU-4NB
Neutralizing Buffer, 5 ml	BLU-5NB
Neutralizing Buffer, 10 ml	BLU-10NB
Dry swab for wet surfaces or add solution	BLU-DRY

Small Tip delivers 1 ml, 100 case	
HiCap™ Neutralizing Broth, 1 ml	BLU-1HC-S
Lethen Broth, 1 ml	BLU-1LET-S

Extensions and Adapters	
PUR-Blue Swab Adapter, 1 unit	BLU-ADAPT
Extension pole, 4 ft, 1 unit	EXT POLE 4FT
Extension pole, 8 ft, 1 unit	EXT POLE 8FT

** Listed in FDA BAM Chapter 10 procedures for *Listeria monocytogenes*



*Add an extension pole
to reach even further*