



Poultry & Carcass Sampling

For Sample Collection



Choose the Right Device

Match the tool to your application



Boot Swabs

Boot swabs are stretchable, sterile fabric boot covers with an integral elastic band designed to easily slide over any shoe or boot and remain securely in place for effortless collection of samples from the floor surface. Available pre-moistened with double-strength skim milk, boot swabs are packaged as singles or pairs in a completely sealed, non-perforated, sterile sample bag to ensure product integrity and quality.

Pull String Sampler

Pull-String samplers are designed for collecting samples from poultry and other livestock housing facilities. The 4" x 4" 12-ply cotton gauze pad with an attached string is available pre-moistened with double-strength skim milk, packaged in a completely sealed, non-perforated, sterile sample bag to ensure product integrity and quality.



Carcass Sampling Kit

The Carcass Sampling Kit contains all components (sponge, gloves, 100 cm² template, and 25 ml of Butterfield's phosphate buffer) needed to collect and transport 15 carcass samples according to the guidelines of the U.S. Pathogen Reduction Act. Kit components (sponges, gloves, template and buffer) can also be bought separately to better match your workflow.



SampleRight™ Sponge Sampler

The SampleRight™ Sponge Sampler with a high-density polyurethane sponge resists crumbling and flaking, making it ideal for carcass sampling. The sponge is easy to manipulate and hold while collecting samples. Available pre-hydrated with Buffered Peptone Water or dry, with or without gloves, to fit your carcass sampling needs.



EZ Reach™ Sponge Sampler

The EZ Reach™ Sponge Sampler with a patented twist-off handle eliminates the need to touch the sponge, reducing the risk of sample contamination and making it a great option for cattle, swine, or even turkey carcass sampling. Available dry, or pre-hydrated with Buffered Peptone water.

The EZ Reach™ 4X sponge sampler, hydrated with double strength skim milk is an excellent alternative for monitoring poultry breeding facilities and grow houses as well as hatcheries and egg laying operations.



US Patent 8,127,627 B2

www.worldbioproducts.com

Collection Devices are just the beginning

Streamline sample collection and processing

Adapters & Extension Poles

Lose the ladder and safely collect samples from hard to reach areas with telescoping extension poles and specialized adapters made exclusively for EZ Reach™ sponges and PUR-Blue™ swabs. Extension poles are available in 4, 8, and 12 ft.



Dividable Grid Containers

Safely and efficiently move samples around the laboratory. Reconfigurable dividers provide clear separation of samples, allowing for easy handling and processing while minimizing the risk of cross contamination.



Online Resource Library

Learn how to improve the effectiveness of your environmental monitoring programs. Visit our online Education Library to download white papers and studies, or view recorded presentations detailing the efficacy and advantages of HiCap Neutralizing Broth.

worldbioproducts.com/education

www.worldbioproducts.com
orders@worldbioproducts.net
PO Box 947, Bothell, WA 98041-0947
(877) 260-6441

WBMKT_2010

