

SampleRight™ Sponge Sampler

Upgrade from traditional cellulose to PUR-polyurethane

Allows for vigorous scrubbing
of surface to disrupt biofilm


Better release of collected microorganisms
with non-toxic, biocide-free polyurethane

Higher density improves rigidity
for easier handling and manipulation

Resists tearing and fraying
unlike cellulose sponges



The SampleRight™ Sponge Sampler has been designed specifically to meet the needs of today's exacting testing procedures. Only SampleRight offers an advance medical grade PUR-polyurethane sponge that resists tearing and flaking, making it the ideal choice for abrasive surfaces and sharp edges. Manufactured to consistent rigorous standards, polyurethane is non-toxic and biocide-free, which greatly improves the recovery and release of collected microorganisms, without any components that may interfere with downstream test methods. Try it with our innovative HiCap™ (High Capacity) Neutralizing Broth, the first collection broth developed specifically for use in the food and beverage industry. By neutralizing residual sanitizers on the sampled surface and maintaining the viability of the collected microorganisms for up to 72 hours, HiCap prevents false negative test results and accurately quantifies microorganisms.



Get Better Results

For More Information

Call: 877-260-6441

Email: info@worldbioproducts.net

| SampleRight™ PUR-polyurethane, 100/case | Item # |
|---|--------------------|
| HiCap™ Neutralizing Broth, 10 ml, 18 oz bag | SR18-10HC-HDPUR |
| HiCap™ Neutralizing Broth, 16 ml, 18 oz bag | SR18-16HC-HDPUR |
| Letheen Broth, 10 ml, 18 oz bag w/gloves | SR18-10LET-HDPUR-G |
| Neutralizing Buffer, 10 ml, 24 oz bag w/gloves | SR-10NB-HDPUR-G |
| Buffered Peptone Water, 20 ml, 18 oz bag w/gloves | SR18-20BPW-HDPUR-G |
| Buffered Peptone Water, 25 ml, 18 oz bag | SR18-25BPW-HDPUR |

SampleRight™ cellulose, 100/case

| | |
|---|--------------|
| HiCap Neutralizing Broth, 10 ml, 18 oz bag | SR18-10HC |
| Letheen Broth, 10 ml, 24 oz bag | SR-10LET |
| Letheen Broth, 10 ml, 24 oz bag w/gloves | SR-10LET-G |
| Letheen Broth, 10 ml, 18 oz bag w/gloves | SR18-10LET-G |
| Neutralizing Buffer, 10 ml, 24 oz bag | SR-10NB |
| Neutralizing Buffer, 10 ml, 24 oz bag w/gloves | SR-10NB-G |
| Neutralizing Buffer, 10 ml, 18oz bag with gloves | SR18-10NB-G |
| Buffered Peptone Water, 10 ml, 24 oz bag w/gloves | SR-10BPW-G |
| Buffered Peptone Water, 10 ml, 18 oz bag | SR18-10BPW |
| Buffered Peptone Water, 10 ml, 18 oz bag w/gloves | SR18-10BPW-G |
| Buffered Peptone Water, 25 ml, 18 oz bag | SR18-25BPW |
| Dry cellulose sponge, 24 oz bag | SR-DRY |
| Dry cellulose sponge, 24 oz bag w/gloves | SR-DRY-G |
| Dry cellulose sponge, 18 oz bag | SR18-DRY |
| Dry cellulose sponge, 18 oz bag w/gloves | SR18-DRY-G |